

Product datasheet for **TA360497**


CHPF Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of human CHPF
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	85 kDa
Gene Name:	chondroitin polymerizing factor
Database Link:	NP_001182660.1 Entrez Gene 79586 Human Q8IZ52
Synonyms:	CHSY2; CSS2; FLJ22678; N-acetylgalactosaminyltransferase
Protein Families:	Druggable Genome
Protein Pathways:	Chondroitin sulfate biosynthesis, Metabolic pathways



[View online »](#)

Product images:

Host: Rabbit
Target Name: CHSS2
Sample Type: Jurkat Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: CHPF
Sample Tissue: Human Jurkat Whole Cell
Antibody Dilution: 1.0 μ g/ml