

## Product datasheet for **TA360500**

### ZNF587 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human ZNF587
Specificity:	<b>Expected reactivity:</b> Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	63 kDa
Gene Name:	zinc finger protein 587
Database Link:	<a href="#">NP_001191746.1</a> <a href="#">Entrez Gene 84914 Human</a> <a href="#">Q96SQ5</a>
Synonyms:	FLJ14710; FLJ20813; FLJ52732; UBF-fl; ZF6
Protein Families:	Transcription Factors



[View online »](#)

**Product images:**


Host: Rabbit  
 Target Name: ZNF587  
 Sample Type: U937 Cell Lysate  
 Antibody Dilution: 1.0µg/ml

Host: Rabbit  
 Target Name: ZNF587  
 Sample Tissue: U937 Whole Cell lysates  
 Antibody Dilution: 1.0ug/ml