

Product datasheet for **TA360512**

RIL (PDLIM4) Rabbit Polyclonal Antibody


Product data:

| | |
|-------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | WB |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | The immunogen is a synthetic peptide directed towards the N terminal region of human PDLIM4 |
| Specificity: | Expected reactivity: Human |
| Formulation: | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i> |
| Concentration: | lot specific |
| Purification: | Affinity purified |
| Conjugation: | Unconjugated |
| Storage: | For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 35 kDa |
| Gene Name: | PDZ and LIM domain 4 |
| Database Link: | NP_001124499.1 Entrez Gene 8572 Human P50479 |
| Background: | This gene encodes a protein which may be involved in bone development. Mutations in this gene are associated with susceptibility to osteoporosis. |
| Synonyms: | RIL |



[View online »](#)

Product images:



Host: Rabbit
Target Name: PDLIM4
Sample Type: Ovary Tumor Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: PDLIM4
Sample Tissue: Human Ovary Tumor
Antibody Dilution: 1.0ug/ml