

Product datasheet for **TA360641**

Zdhhc4 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C terminal region of rat ZDHHC4
Specificity:	Expected reactivity: Rat
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	39 kDa
Gene Name:	zinc finger, DHHC-type containing 4
Database Link:	NP_001013141.1 Entrez Gene 304291 Rat F1LQN6
Synonyms:	DHHC-4; FLJ10479; OTTHUMP00000115925; OTTHUMP00000115926; ZNF374



[View online »](#)

Product images:



Host: Rabbit
Target Name: Zdhhc4
Sample Type: Rat Muscle Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: ZDHHC4
Sample Tissue: Rat Muscle lysates
Antibody Dilution: 1.0ug/ml