

Product datasheet for **TA360692**

GLIPR2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human GLIPR2
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	16 kDa
Gene Name:	GLI pathogenesis related 2
Database Link:	NP_071738.1 Entrez Gene 152007 Human Q9H4G4
Synonyms:	C9orf19; GAPR-1; GAPR1



[View online »](#)

Product images:



Host: Rabbit
Target Name: GLIPR2
Sample Type: MCF7 Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: GLIPR2
Sample Tissue: MCF7 Whole Cell lysates
Antibody Dilution: 1.0ug/ml