

Product datasheet for **TA360708**

OSBPL10 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of human OSBPL10
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	84 kDa
Gene Name:	oxysterol binding protein like 10
Database Link:	NP_001167531.1 Entrez Gene 114884 Human Q9BXB5
Background:	This gene encodes a member of the oxysterol-binding protein (OSBP) family, a group of intracellular lipid receptors. Like most members, the encoded protein contains an N-terminal pleckstrin homology domain and a highly conserved C-terminal OSBP-like sterol-binding domain. Multiple transcript variants encoding different isoforms have been found for this gene.
Synonyms:	FLJ20363; ORP-10; ORP10; OSBP9; OTTHUMP00000161257



[View online »](#)

Product images:



Host: Rabbit
Target Name: OSBPL10
Sample Type: Thyroid Tumor Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: OSBPL10
Sample Tissue: Human Thyroid Tumor lysates
Antibody Dilution: 1 ug/ml