

Product datasheet for **TA360728**

MFSD6 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of human MFSD6
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	87 kDa
Gene Name:	major facilitator superfamily domain containing 6
Database Link:	NP_060164.3 Entrez Gene 54842 Human Q6ZSS7
Synonyms:	FLJ20160; MMR2; OTTHUMP00000163510
Protein Families:	Transmembrane



[View online »](#)

Product images:



Host: Rabbit
Target Name: MFSD6
Sample Type: A172 Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: MFSD6
Sample Tissue: A172 Whole Cell lysates
Antibody Dilution: 1.0ug/ml