

## Product datasheet for **TA360755**


### SPATA25 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human SPATA25
Specificity:	<b>Expected reactivity:</b> Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	24 kDa
Gene Name:	spermatogenesis associated 25
Database Link:	<a href="#">NP_542175.1</a> <a href="#">Entrez Gene 128497 Human</a> <a href="#">Q9BR10</a>
Synonyms:	1700020C07Rik; C20orf165; RP23-61O3.3; Tsg23

[View online »](#)

## Product images:



Host: Rabbit  
Target Name: C20orf165  
Sample Type: Lung Tumor Lysate  
Antibody Dilution: 1.0µg/ml

Host: Rabbit  
Target Name: C20ORF165  
Sample Tissue: Lung Tumor lysates  
Antibody Dilution: 1ug/ml