

Product datasheet for TA360777

SLC47A1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human SLC47A1
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2–8°C up to 1 week. For long term storage, store at –20°C in small aliquots to prevent freeze–thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	29 kDa
Gene Name:	solute carrier family 47 member 1
Database Link:	NP_060712.2 Entrez Gene 55244 Human I3L345
Background:	This gene is located within the Smith–Magenis syndrome region on chromosome 17. It encodes a protein of unknown function.



Synonyms: FLJ10847; hMATE-1; MATE-1; MATE1; MGC64822

Product images:



Host: Rabbit
Target Name: SLC47A1
Sample Type: 786-0 Cell Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: SLC47A1
Sample Tissue: Human 786-0 Whole Cell
lysates
Antibody Dilution: 1ug/ml