

Product datasheet for TA360840

SH3BGRL3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human SH3BGRL3
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2–8°C up to 1 week. For long term storage, store at –20°C in small aliquots to prevent freeze–thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	10 kDa
Gene Name:	SH3 domain binding glutamate rich protein like 3
Database Link:	NP_112576.1 Entrez Gene 83442 Human Q9H299
Background:	Could act as a modulator of glutaredoxin biological activity.



Synonyms: SH3BP-1; TIP-B1

Product images:



Host: Rabbit
Target Name: SH3BGRL3
Sample Type: Ovary Tumor Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: SH3BGRL3
Sample Tissue: Human Ovary Tumor lysates
Antibody Dilution: 1ug/ml