

Product datasheet for **TA360863**

PAN3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human PAN3
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	63 kDa
Gene Name:	PAN3 poly(A) specific ribonuclease subunit
Database Link:	NP_787050.6 Entrez Gene 255967 Human Q58A45-2
Synonyms:	hPan3; OTTHUMP00000018179
Protein Families:	Druggable Genome



[View online »](#)

Product images:



Host: Rabbit
Target Name: PAN3
Sample Type: Stomach Tumor Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: PAN3
Sample Tissue: Human Stomach Tumor lysates
Antibody Dilution: 1 μ g/ml