

Product datasheet for **TA360900**

TDIF2 (DNTTIP2) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human DNTTIP2
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	83 kDa
Gene Name:	deoxynucleotidyltransferase terminal interacting protein 2
Database Link:	NP_055412.2 Entrez Gene 30836 Human P03372
Background:	This gene is thought to be involved in chromatin remodeling and gene transcription. The encoded nuclear protein binds to and enhances the transcriptional activity of the estrogen receptor alpha, and also interacts with terminal deoxynucleotidyltransferase. The expression profile of this gene is a potential biomarker for chronic obstructive pulmonary disease.
Synonyms:	ERBP; FCF2; HSU15552; LPTS-RP2; MGC163494; RP4-561L24.1; TdIF2
Protein Families:	Druggable Genome



[View online »](#)

Product images:



Host: Rabbit
Target Name: DNTTIP2
Sample Type: COLO205 Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: DNTTIP2
Sample Tissue: Human COLO205 Whole Cell lysates
Antibody Dilution: 1 μ g/ml