

Product datasheet for TA360902


FAM116A (DENND6A) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human FAM116A
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	70 kDa
Gene Name:	DENN domain containing 6A
Database Link:	NP_689891.1 Entrez Gene 201627 Human Q8IWF6
Synonyms:	AF11A; FAM116A



Product images:



Host: Rabbit
 Target Name: FAM116A
 Sample Type: 721_B Cell Lysate
 Antibody Dilution: 1.0µg/ml

Host: Rabbit
 Target Name: FAM116A
 Sample Tissue: Human 721_B lymphoblast
 Whole Cell lysates
 Antibody Dilution: 1ug/ml