

Product datasheet for **TA360927**

MOBKL2A (MOB3A) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of human MOBKL2A
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	24 kDa
Gene Name:	MOB kinase activator 3A
Database Link:	NP_570719.1 Entrez Gene 126308 Human Q96BX8
Synonyms:	5330417K06Rik; A630029F06; AV218468; Mobkl2a; Mobkl2b



[View online »](#)

Product images:



Host: Rabbit
Target Name: MOBKL2A
Sample Type: U937 Cell Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: MOBKL2A
Sample Tissue: Human U937 Whole Cell lysates
Antibody Dilution: 1ug/ml