

Product datasheet for **TA360959**

Huntingtin Associated Protein 1 (HAP1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human HAP1
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	74 kDa
Gene Name:	huntingtin-associated protein 1
Database Link:	NP_001073339.1 Entrez Gene 9001 Human P54257



[View online »](#)

Background:


Huntington's disease (HD), a neurodegenerative disorder characterized by loss of striatal neurons, is caused by an expansion of a polyglutamine tract in the HD protein huntingtin. This gene encodes a protein that interacts with huntingtin, with two cytoskeletal proteins (dynactin and pericentriolar autoantigen protein 1), and with a hepatocyte growth factor-regulated tyrosine kinase substrate. The interactions with cytoskeletal proteins and a kinase substrate suggest a role for this protein in vesicular trafficking or organelle transport. Several alternatively spliced transcript variants encoding different isoforms have been described for this gene.

Synonyms:

HAP-1; HAP2; hHLP1; HIP5; HLP; HLP1

Protein Pathways:

Huntington's disease

Product images:


Host: Rabbit
 Target Name: HAP1
 Sample Type: U937 Cell Lysate
 Antibody Dilution: 1.0µg/ml

Host: Rabbit
 Target Name: HAP1
 Sample Tissue: Human U937 Whole Cell lysates
 Antibody Dilution: 1ug/ml