

Product datasheet for **TA360961**

HCCA2 (MOB2) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human MOB2
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	26 kDa
Gene Name:	MOB kinase activator 2
Database Link:	NP_001165694.1 Entrez Gene 81532 Human Q70IA6
Background:	Stimulates the autophosphorylation and kinase activity of STK38 and STK38L.
Synonyms:	HCCA2



[View online »](#)

Product images:



Host: Rabbit
Target Name: MOB2
Sample Type: PANC1 Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: MOB2
Sample Tissue: Human PANC1 Whole Cell lysates
Antibody Dilution: 1 μ g/ml