

## Product datasheet for **TA360985**

### FBRS Rabbit Polyclonal Antibody


#### Product data:

|                                |                                                                                                                                                                            |
|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b>           | Primary Antibodies                                                                                                                                                         |
| <b>Applications:</b>           | WB                                                                                                                                                                         |
| <b>Reactivity:</b>             | Human                                                                                                                                                                      |
| <b>Host:</b>                   | Rabbit                                                                                                                                                                     |
| <b>Clonality:</b>              | Polyclonal                                                                                                                                                                 |
| <b>Immunogen:</b>              | The immunogen is a synthetic peptide directed towards the middle region of human FBRS                                                                                      |
| <b>Specificity:</b>            | <b>Expected reactivity:</b> Human                                                                                                                                          |
| <b>Formulation:</b>            | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.                                                                          |
| <b>Concentration:</b>          | lot specific                                                                                                                                                               |
| <b>Purification:</b>           | Affinity purified                                                                                                                                                          |
| <b>Conjugation:</b>            | Unconjugated                                                                                                                                                               |
| <b>Storage:</b>                | For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.                                    |
| <b>Stability:</b>              | Shelf life: one year from despatch.                                                                                                                                        |
| <b>Predicted Protein Size:</b> | 41 kDa                                                                                                                                                                     |
| <b>Gene Name:</b>              | fibrosin                                                                                                                                                                   |
| <b>Database Link:</b>          | <a href="#">NP_001098549.2</a><br><a href="#">Entrez Gene 64319 Human</a><br><a href="#">Q9HAH7</a>                                                                        |
| <b>Background:</b>             | Fibrosin is a lymphokine secreted by activated lymphocytes that induces fibroblast proliferation (Prakash and Robbins, 1998 [PubMed 9809749]).[supplied by OMIM, Mar 2008] |
| <b>Synonyms:</b>               | FBS; FBS1; fibrosin; FLJ11618                                                                                                                                              |



Protein Families: Secreted Protein

### Product images:



Host: Rabbit  
Target Name: FBR5  
Sample Type: Jurkat Cell Lysate  
Antibody Dilution: 1.0µg/ml

Host: Rabbit  
Target Name: FBR5  
Sample Tissue: Human Jurkat Whole Cell  
lysates  
Antibody Dilution: 1ug/ml