

Product datasheet for **TA361030**

ACSM2B Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C terminal region of human ACSM2B
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	63 kDa
Gene Name:	acyl-CoA synthetase medium-chain family member 2B
Database Link:	NP_001098539.1 Entrez Gene 348158 Human Q68CK6
Synonyms:	ACSM2; HXMA; HYST1046



[View online »](#)

Product images:



Host: Rabbit
Target Name: ACSM2B
Sample Type: Lung Tumor Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: ACSM2B
Sample Tissue: Human Lung Tumor lysates
Antibody Dilution: 1ug/ml