

Product datasheet for TA361147

BPIFC Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human BPIL2
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	55 kDa
Gene Name:	BPI fold containing family C
Database Link:	NP_777592.1 Entrez Gene 254240 Human Q8NFAQ6
Synonyms:	BPIL2
Protein Families:	Druggable Genome



Product images:



Host: Rabbit
Target Name: BPIL2
Sample Type: HeLa Cell Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: BPIL2
Sample Tissue: Human HeLa Whole Cell
lysates
Antibody Dilution: 1ug/ml