

Product datasheet for **TA361317**

SMC6L1 (SMC6) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of human SMC6
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	120 kDa
Gene Name:	structural maintenance of chromosomes 6
Database Link:	NP_001135758.1 Entrez Gene 79677 Human Q96SB8




[View online »](#)

Background:

Core component of the SMC5-SMC6 complex, a complex involved in DNA double-strand breaks by homologous recombination. The complex may promote sister chromatid homologous recombination by recruiting the SMC1-SMC3 cohesin complex to double-strand breaks. The complex is required for telomere maintenance via recombination in ALT (alternative lengthening of telomeres) cell lines and mediates sumoylation of shelterin complex (telosome) components which is proposed to lead to shelterin complex disassembly in ALT-associated PML bodies (APBs). Required for recruitment of telomeres to PML nuclear bodies. SMC5-SMC6 complex may prevent transcription of episomal DNA, such as circular viral DNA genome.

Synonyms:

FLJ22116; FLJ35534; hSMC6; OTTHUMP00000200490; SMC6L1

Product images:

Host: Rabbit
Target Name: SMC6
Sample Type: Neurofibroma Tumor Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: SMC6
Sample Tissue: Human Neurofibroma Tumor lysates
Antibody Dilution: 1ug/ml