

## Product datasheet for **TA361547**

### RPH3AL Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human RPH3AL
Specificity:	<b>Expected reactivity:</b> Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	34 kDa
Gene Name:	rabphilin 3A-like (without C2 domains)
Database Link:	<a href="#">NP_001177340.1</a> <a href="#">Entrez Gene 9501 Human</a> <a href="#">Q9UNE2</a>
Background:	The protein encoded by this gene plays a direct regulatory role in calcium-ion-dependent exocytosis in both endocrine and exocrine cells and plays a key role in insulin secretion by pancreatic cells. This gene is likely a tumor suppressor. Alternative splicing results in multiple transcript variants encoding distinct isoforms.
Synonyms:	NOC2
Protein Families:	Secreted Protein



[View online »](#)

## Product images:



Host: Rabbit

Target Name: RPH3AL

Sample Type: Leiomyosarcoma Tumor Lysate

Antibody Dilution: 1.0 $\mu$ g/ml

Host: Rabbit

Target Name: RPH3AL

Sample Tissue: Human Leiomyosarcoma Tumor lysates

Antibody Dilution: 1 $\mu$ g/ml