

Product datasheet for **TA361580**

DNAJA4 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of human DNAJA4
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	43 kDa
Gene Name:	DnaJ heat shock protein family (Hsp40) member A4
Database Link:	NP_001123654.1 Entrez Gene 55466 Human Q8WW22
Background:	The function of this protein remains unknown.
Synonyms:	MST104; MSTP104; PRO1472



[View online »](#)

Product images:



Host: Rabbit
Target Name: DNAJA4
Sample Type: HCT116 Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: DNAJA4
Sample Tissue: Human HCT116 Whole Cell lysates
Antibody Dilution: 1 μ g/ml