

## Product datasheet for **TA361594**

### TNS4 Rabbit Polyclonal Antibody


#### Product data:

|                         |   |
|-------------------------|---|
| Product Type:           | Primary Antibodies  |
| Applications:           | WB  |
| Reactivity:             | Human   |
| Host:                   | Rabbit  |
| Clonality:              | Polyclonal  |
| Immunogen:              | The immunogen is a synthetic peptide directed towards the middle region of human TNS4   |
| Specificity:            | <b>Expected reactivity:</b> Human   |
| Formulation:            | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.<br><i>Note that this product is shipped as lyophilized powder to China customers.</i> |
| Concentration:          | lot specific  |
| Purification:           | Affinity purified   |
| Conjugation:            | Unconjugated  |
| Storage:                | For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.   |
| Stability:              | Shelf life: one year from despatch.   |
| Predicted Protein Size: | 78 kDa  |
| Gene Name:              | tensin 4  |
| Database Link:          | <a href="#">NP_116254.4</a><br><a href="#">Entrez Gene 84951 Human</a><br><a href="#">Q8IZW8</a>  |
| Background:             | May be involved in cell migration, cartilage development and in linking signal transduction pathways to the cytoskeleton. May promote apoptosis, via its cleavage by caspase-3.         |
| Synonyms:               | CTEN; FLJ14950; OTTHUMP00000164465; PP14434   |



[View online »](#)

## Product images:



Host: Rabbit  
Target Name: TNS4  
Sample Type: HeLa Cell Lysate  
Antibody Dilution: 1.0 $\mu$ g/ml

Host: Rabbit  
Target Name: TNS4  
Sample Tissue: Human HeLa Whole Cell lysates  
Antibody Dilution: 1 $\mu$ g/ml