

Product datasheet for **TA361611**


CTRP1 (C1QTNF1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human C1QTNF1
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	30 kDa
Gene Name:	C1q and tumor necrosis factor related protein 1
Database Link:	NP_112230.1 Entrez Gene 114897 Human Q9BXJ1
Background:	The function of this protein remains unknown.
Synonyms:	CTRP1; FLJ90694; GIP; ZSIG37
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane



[View online »](#)

Product images:

Host: Rabbit
Target Name: C1QTNF1
Sample Type: Testicular Tumor Lysate
Antibody Dilution: 1.0 µg/ml

Host: Rabbit
Target Name: C1QTNF1
Sample Tissue: Human Testicular Tumor lysates
Antibody Dilution: 1 µg/ml