

## Product datasheet for **TA361714**

### Tetraspanin 9 (TSPAN9) Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C terminal region of human TSPAN9
Specificity:	<b>Expected reactivity:</b> Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	26 kDa
Gene Name:	tetraspanin 9
Database Link:	<a href="#">NP_001161792.1</a> <a href="#">Entrez Gene 10867 Human</a> <a href="#">Q75954</a>
Background:	The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. Alternatively spliced transcripts encoding the same protein have been identified.
Synonyms:	NET-5; NET5; PP1057; Tspan-9
Protein Families:	Transmembrane



[View online »](#)

## Product images:



Host: Rabbit  
Target Name: TSPAN9  
Sample Type: HepG2 Cell Lysate  
Antibody Dilution: 1.0µg/ml

Host: Rabbit  
Target Name: TSPAN9  
Sample Tissue: Human HepG2 Whole Cell lysates  
Antibody Dilution: 1ug/ml