

Product datasheet for **TA361759**

Methyltransferase like 5 (METTL5) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle terminal region of human METTL5
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	22 kDa
Gene Name:	methyltransferase like 5
Database Link:	NP_001280115.1 Entrez Gene 29081 Human Q9NRN9
Background:	Probable methyltransferase.
Synonyms:	FLJ10459; HSPC133; OTTHUMP00000205106
Protein Families:	Druggable Genome



[View online »](#)

Product images:



Host: Rabbit
Target Name: METTL5
Sample Type: Ovary Tumor Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: METTL5
Sample Tissue: Human Ovary Tumor lysates
Antibody Dilution: 1ug/ml