

Product datasheet for **TA361777**

DYNC112 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human DYNC112
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	70 kDa
Gene Name:	dynein cytoplasmic 1 intermediate chain 2
Database Link:	NP_001258714.1 Entrez Gene 1781 Human Q13409
Background:	This gene encodes a member of the dynein intermediate chain family. The encoded protein is a non-catalytic component of the cytoplasmic dynein 1 complex, which acts as a retrograde microtubule motor to transport organelles and vesicles. A pseudogene of this gene is located on chromosome 10. Alternative splicing results in multiple transcript variants.
Synonyms:	DNCI2; DNCIC2; FLJ21089; FLJ90842; IC2; MGC9324; MGC104199; MGC111094; OTTHUMP00000205409



[View online »](#)

Product images:



Host: Rabbit
Target Name: DYNC112
Sample Type: Stomach Tumor Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: DYNC112
Sample Tissue: Human Stomach Tumor lysates
Antibody Dilution: 1 μ g/ml