

Product datasheet for **TA362085**

BRCC36 (BRCC3) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle terminal region of human BRCC3
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	36 kDa
Gene Name:	BRCA1/BRCA2-containing complex subunit 3
Database Link:	NP_001018065.1 Entrez Gene 79184 Human P46736



[View online »](#)

Background:


This gene encodes a subunit of the BRCA1-BRCA2-containing complex (BRCC), which is an E3 ubiquitin ligase. This complex plays a role in the DNA damage response, where it is responsible for the stable accumulation of BRCA1 at DNA break sites. The component encoded by this gene can specifically cleave Lys 63-linked polyubiquitin chains, and it regulates the abundance of these polyubiquitin chains in chromatin. The loss of this gene results in abnormal angiogenesis and is associated with syndromic moyamoya, a cerebrovascular angiopathy. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 5.

Synonyms:

BRCC36; c6.1A; CXorf53

Protein Families:

Druggable Genome, Protease

Product images:


Host: Rabbit
 Target Name: BRCC3
 Sample Type: Du145 Cell Lysate
 Antibody Dilution: 1.0µg/ml

Host: Rabbit
 Target Name: BRCC3
 Sample Tissue: Human Du145 Whole Cell lysates
 Antibody Dilution: 1ug/ml