

Product datasheet for **TA362227**


Centaurin beta 2 (ACAP2) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human ACAP2
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	85 kDa
Gene Name:	ArfGAP with coiled-coil, ankyrin repeat and PH domains 2
Database Link:	NP_036419.3 Entrez Gene 23527 Human Q15057
Background:	GTPase-activating protein (GAP) for ADP ribosylation factor 6 (ARF6).
Synonyms:	Centaurin-beta-2; CENTB2; CNT-B2; KIAA0041
Protein Pathways:	Endocytosis



[View online »](#)

Product images:

Host: Rabbit
Target Name: ACAP2
Sample Type: Jurkat Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: ACAP2
Sample Tissue: Human Jurkat Whole Cell lysates
Antibody Dilution: 1 μ g/ml