

## Product datasheet for **TA362268**

### DHX15 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of human DHX15
Specificity:	<b>Expected reactivity:</b> Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	87 kDa
Gene Name:	DEAH-box helicase 15
Database Link:	<a href="#">NP_001349.2</a> <a href="#">Entrez Gene 1665 Human</a> <a href="#">O43143</a>
Background:	The protein encoded by this gene is a putative ATP-dependent RNA helicase implicated in pre-mRNA splicing.
Synonyms:	DBP1; DDX15; HRH2; PRP43; PRPF43; PrPp43p
Protein Pathways:	Spliceosome



[View online »](#)

## Product images:



Host: Rabbit  
Target Name: DHX15  
Sample Type: 786-0 Cell Lysate  
Antibody Dilution: 1.0 $\mu$ g/ml

Host: Rabbit  
Target Name: DHX15  
Sample Tissue: Human 786-0 Whole Cell lysates  
Antibody Dilution: 1 $\mu$ g/ml