

Product datasheet for **TA362285**

RAB36 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human RAB36
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	36 kDa
Gene Name:	RAB36, member RAS oncogene family
Database Link:	NP_004905.2 Entrez Gene 9609 Human O95755
Background:	Protein transport. Probably involved in vesicular traffic (By similarity).
Synonyms:	RAB36
Protein Families:	Druggable Genome



[View online »](#)

Product images:



Host: Rabbit
Target Name: RAB36
Sample Type: Du145 Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: RAB36
Sample Tissue: Human Du145 Whole Cell lysates
Antibody Dilution: 1 ug/ml