

Product datasheet for TA362428


THAP4 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human THAP4
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2–8°C up to 1 week. For long term storage, store at –20°C in small aliquots to prevent freeze–thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	63 kDa
Gene Name:	THAP domain containing 4
Database Link:	NP_001157828.1 Entrez Gene 51078 Human Q8WY91
Background:	The function of this protein remains unknown.
Synonyms:	CGI-36; PP238



Product images:



Host: Rabbit
Target Name: THAP4
Sample Type: RPMI 8226 Cell Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: THAP4
Sample Tissue: Human RPMI 8226 Whole Cell
lysates
Antibody Dilution: 1ug/ml