

Product datasheet for **TA362461**


LRRC40 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C terminal region of human LRRC40
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	66 kDa
Gene Name:	leucine rich repeat containing 40
Database Link:	NP_060238.3 Entrez Gene 55631 Human Q9H9A6
Background:	The function of this protein remains unknown.
Synonyms:	dj677H15.1; FLJ20331; RP4-677H15.1
Protein Families:	Druggable Genome



[View online »](#)

Product images:

Host: Rabbit
Target Name: LRRC40
Sample Type: HeLa Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: LRRC40
Sample Tissue: Human HeLa Whole Cell lysates
Antibody Dilution: 1 μ g/ml