

Product datasheet for **TA362475**

MINDY1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C terminal region of human FAM63A
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	57 kDa
Gene Name:	family with sequence similarity 63 member A
Database Link:	NP_001035307.1 Entrez Gene 55793 Human Q8N5J2-3
Synonyms:	FLJ11280; FLJ43504; KIAA1390; RP11-316M1.5



[View online »](#)

Product images:



Host: Rabbit
Target Name: FAM63A
Sample Type: Jurkat Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: FAM63A
Sample Tissue: Human Jurkat Whole Cell lysates
Antibody Dilution: 1 μ g/ml