

Product datasheet for **TA362487**

UBXN6 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human UBXN6
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	48 kDa
Gene Name:	UBX domain protein 6
Database Link:	NP_001164562.1 Entrez Gene 80700 Human Q9BZV1
Background:	Acts in a complex with VCP and cooperates with USP7 in promoting MDM2 deubiquitination and stabilization. MDM2 stabilization leads to MDM2-dependent TP53 degradation.
Synonyms:	DKFZp667D109; UBXD1; UBXDC2



[View online »](#)

Product images:



Host: Rabbit
 Target Name: UBXN6
 Sample Type: Jurkat Cell Lysate
 Antibody Dilution: 1.0µg/ml

Host: Rabbit
 Target Name: UBXN6
 Sample Tissue: Human Jurkat Whole Cell lysates
 Antibody Dilution: 1ug/ml