

Product datasheet for **TA362666**

MYRIP Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle terminal region of human MYRIP
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	87 kDa
Gene Name:	myosin VIIA and Rab interacting protein
Database Link:	NP_001271352.1 Entrez Gene 25924 Human Q8NFW9-2
Background:	Rab effector protein involved in melanosome transport. Serves as link between melanosome-bound RAB27A and the motor proteins MYO5A and MYO7A. May link RAB27A-containing vesicles to actin filaments. Functions as a protein kinase A-anchoring protein (AKAP). May act as a scaffolding protein that links PKA to components of the exocytosis machinery, thus facilitating exocytosis, including insulin release (By similarity).
Synonyms:	DKFZp586F1018; exophilin-8; FLJ44025; MGC130034; MGC130035; OTTHUMP00000209206; SLAC2-C; SLAC2C



[View online »](#)

Product images:



Host: Rabbit
Target Name: MYRIP
Sample Type: PANC1 Cell Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: MYRIP
Sample Tissue: Human PANC1 Whole Cell lysates
Antibody Dilution: 1ug/ml