

Product datasheet for TA362725

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

HIP1 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

WB **Applications:**

Reactivity: Human Rabbit Host:

Clonality: Polyclonal

Immunogen: The immunogen is a synthetic peptide directed towards the middle region of human HIP1

Specificity: **Expected reactivity**: Human

Formulation: Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2%

Note that this product is shipped as lyophilized powder to China customers.

Concentration: lot specific

Purification: Affinity purified Conjugation: Unconjugated

For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small Storage:

aliquots to prevent freeze-thaw cycles.

Stability: Shelf life: one year from despatch.

Predicted Protein Size: 116 kDa

Gene Name: huntingtin interacting protein 1

Database Link: NP 001230127.1

Entrez Gene 3092 Human

O00291

Background: The product of this gene is a membrane-associated protein that functions in clathrin-

> mediated endocytosis and protein trafficking within the cell. The encoded protein binds to the huntingtin protein in the brain; this interaction is lost in Huntington's disease. Alternative

splicing results in multiple transcript variants.

HIP-I; ILWEQ; MGC126506 Synonyms:


Protein Families: Druggable Genome





Protein Pathways: Huntington's disease

Product images:



Host: Rabbit Target Name: HIP1

Sample Tissue: Human HepG2 Whole Cell lysates

Antibody Dilution: 1ug/ml