

Product datasheet for **TA362773**

Ku70 (XRCC6) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C terminal region of human XRCC6
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	66 kDa
Gene Name:	X-ray repair complementing defective repair in Chinese hamster cells 6
Database Link:	NP_001275905.1 Entrez Gene 2547 Human P12956
Background:	The p70/p80 autoantigen is a nuclear complex consisting of two subunits with molecular masses of approximately 70 and 80 kDa. The complex functions as a single-stranded DNA-dependent ATP-dependent helicase. The complex may be involved in the repair of nonhomologous DNA ends such as that required for double-strand break repair, transposition, and V(D)J recombination. High levels of autoantibodies to p70 and p80 have been found in some patients with systemic lupus erythematosus.



[View online »](#)

Synonyms: CTC75; CTCBF; D22S671; D22S731; G22P1; KU70; ML8; OTTHUMP00000199333; OTTHUMP00000199349; TLA

Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Non-homologous end-joining

Product images:



Host: Rabbit

Target Name: XRCC6

Sample Type: Jurkat Cell Lysate

Antibody Dilution: 1.0 μ g/ml

Host: Rabbit

Target Name: XRCC6

Sample Tissue: Human Jurkat Whole Cell lysates

Antibody Dilution: 1 μ g/ml