

Product datasheet for **TA362846**


USP33 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human USP33
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	103 kDa
Gene Name:	ubiquitin specific peptidase 33
Database Link:	NP_055832.3 Entrez Gene 23032 Human Q8TEY7
Background:	This gene encodes a deubiquinating enzyme important in a variety of processes, including Slit-dependent cell migration and beta-2 adrenergic receptor signaling. The protein is negatively regulated through ubiquitination by von Hippel-Lindau tumor protein (VHL). Alternative splicing results in multiple transcript variants and protein isoforms.
Synonyms:	hVDU1; KIAA1097; MGC16868; OTTHUMP00000011330; VDU1
Protein Families:	Druggable Genome, Protease



[View online »](#)

Product images:

Host: Rabbit
Target Name: USP33
Sample Type: THP-1 Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: USP33
Sample Tissue: Human THP-1 Whole Cell lysates
Antibody Dilution: 1 μ g/ml