

Product datasheet for **TA362943**

Ugt1a1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of mouse UGT1A1
Specificity:	Expected reactivity: Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	58 kDa
Gene Name:	UDP glucuronosyltransferase 1 family, polypeptide A1
Database Link:	NP_964007.2 Entrez Gene 394436 Mouse Q63886
Background:	UDPGT is of major importance in the conjugation and subsequent elimination of potentially toxic xenobiotics and endogenous compounds.
Synonyms:	GNT1; HUG-BR1; UDPGT; UGT-1A; UGT1; UGT1*1; UGT1-01; UGT1.1; UGT1A



[View online »](#)

Product images:



Host: Rabbit
Target Name: Ugt1a1
Sample Type: Mouse Kidney Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: UGT1A1
Sample Tissue: Mouse Kidney lysates
Antibody Dilution: 1 μ g/ml