

## Product datasheet for **TA362963**

### Pecam1 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of mouse PECAM1
Specificity:	<b>Expected reactivity:</b> Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	68 kDa
Gene Name:	platelet and endothelial cell adhesion molecule 1
Database Link:	<a href="#">NP_113779.1</a> <a href="#">Q08481-4</a>
Background:	cell adhesion molecule, also intracellular signaling molecule; may mediate leukocyte adhesion and migration, angiogenesis, and thrombosis.
Synonyms:	CD31; CD31/EndoCAM; EndoCAM; FLJ58394; PECAM-1



[View online »](#)

## Product images:



Host: Rabbit  
Target Name: Pecam1  
Sample Type: Mouse Brain Lysate  
Antibody Dilution: 1.0 $\mu$ g/ml

Host: Rabbit  
Target Name: PECAM1  
Sample Tissue: Mouse Brain lysates  
Antibody Dilution: 1 $\mu$ g/ml