

Product datasheet for **TA363049**

Pdia3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of mouse PDIA3
Specificity:	Expected reactivity: Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	55 kDa
Gene Name:	protein disulfide isomerase associated 3
Database Link:	NP_031978.2 Entrez Gene 14827 Mouse P27773
Background:	The function of this protein remains unknown.
Synonyms:	ER60; ERp57; ERp60; ERp61; GRP57; GRP58; HsT17083; P58; PI-PLC



[View online »](#)

Product images:



Host: Rabbit
Target Name: Pdia3
Sample Type: Mouse Heart Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: PDIA3
Sample Tissue: Mouse Heart lysates
Antibody Dilution: 1ug/ml