

Product datasheet for **TA363052**

Fabp3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of mouse FABP3
Specificity:	Expected reactivity: Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	14 kDa
Gene Name:	fatty acid binding protein 3, muscle and heart
Database Link:	NP_034304.1 Entrez Gene 14077 Mouse P11404
Background:	FABP are thought to play a role in the intracellular transport of long-chain fatty acids and their acyl-CoA esters.
Synonyms:	FABP11; H-FABP; M-FABP; MDGI; O-FABP



[View online »](#)

Product images:



Host: Rabbit
Target Name: Fabp3
Sample Type: Mouse Kidney Lysate
Antibody Dilution: 1.0 µg/ml

Host: Rabbit
Target Name: FABP3
Sample Tissue: Mouse Kidney lysates
Antibody Dilution: 1 µg/ml