

Product datasheet for **TA363102**

Gsta4 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of mouse GSTA4
Specificity:	Expected reactivity: Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	24 kDa
Gene Name:	glutathione S-transferase, alpha 4
Database Link:	NP_034487.2 Entrez Gene 14860 Mouse P24472
Background:	Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.
Synonyms:	DKFZp686D21185; GSTA4-4; GTA4; OTTHUMP00000016625



[View online »](#)

Product images:



Host: Rabbit
Target Name: Gsta4
Sample Type: Mouse Testis Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: GSTA4
Sample Tissue: Mouse Testis lysates
Antibody Dilution: 1 μ g/ml