

Product datasheet for **TA363130**

Acadl Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C terminal region of mouse ACADL
Specificity:	Expected reactivity: Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	47 kDa
Gene Name:	acyl-Coenzyme A dehydrogenase, long-chain
Database Link:	NP_031407.2 Entrez Gene 11363 Mouse P51174




[View online »](#)

Background:

This gene encodes a homotetrameric mitochondrial flavoprotein and is a member of the acyl-CoA dehydrogenase family. Members of this family catalyze the first step of fatty acid beta-oxidation, forming a C2-C3 trans-double bond in a FAD-dependent reaction. As beta-oxidation cycles through its four steps, each member of the acyl-CoA dehydrogenase family works at an optimum fatty acid chain-length. This enzyme has its optimum length between C12- and C16-acylCoA. In mice, deficiency of this gene can cause sudden death, cardiomyopathy as well as fasting and cold intolerance.

Synonyms:

ACAD4; FLJ94052; LCAD

Product images:

Host: Rabbit
Target Name: Acadl
Sample Type: Mouse Small Intestine Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: ACADL
Sample Tissue: Mouse Small Intestine lysates
Antibody Dilution: 1 ug/ml