

Product datasheet for **TA363200**

Klf3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of mouse KLF3
Specificity:	Expected reactivity: Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	37 kDa
Gene Name:	Kruppel-like factor 3 (basic)
Database Link:	NP_032479.1 Entrez Gene 16599 Mouse Q60980
Background:	Binds to the CACCC box of erythroid cell-expressed genes. May play a role in hematopoiesis.
Synonyms:	BKLF; MGC48279; TEF-2



[View online »](#)

Product images:



Host: Rabbit
Target Name: Klf3
Sample Type: Mouse Lung Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: KLF3
Sample Tissue: Mouse Lung lysates
Antibody Dilution: 1 μ g/ml