

## Product datasheet for **TA363207**

### Acaa2 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of mouse ACAA2
Specificity:	<b>Expected reactivity:</b> Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	42 kDa
Gene Name:	acetyl-Coenzyme A acyltransferase 2 (mitochondrial 3-oxoacyl-Coenzyme A thiolase)
Database Link:	<a href="#">NP_803421.1</a> <a href="#">Entrez Gene 52538 Mouse</a> <a href="#">Q8BWT1</a>
Background:	Abolishes BNIP3-mediated apoptosis and mitochondrial damage.
Synonyms:	Beta-ketothiolase; DSAEC; FLJ35992; FLJ95265; T1



[View online »](#)

## Product images:



Host: Rabbit  
Target Name: Acaa2  
Sample Type: Mouse Testis Lysate  
Antibody Dilution: 1.0 $\mu$ g/ml

Host: Rabbit  
Target Name: ACAA2  
Sample Tissue: Mouse Testis lysates  
Antibody Dilution: 1 $\mu$ g/ml