

Product datasheet for **TA363299**

Siah1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C terminal region of rat SIAH1A
Specificity:	Expected reactivity: Rat
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	31 kDa
Gene Name:	siah E3 ubiquitin protein ligase 1
Database Link:	NP_543181.2 Entrez Gene 140941 Rat Q920M9
Background:	synaptophysin binding protein and function as E3 ubiquitin-protein ligases to regulate the ubiquitination and degradation of synaptophysin [RGD, Feb 2006]
Synonyms:	FLJ08065; hSIAH1; HUMSIAH; Siah-1; Siah-1a; SIAH1A



[View online »](#)

Product images:



Host: Rabbit
Target Name: Siah1a
Sample Type: Rat Spleen Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: SIAH1A
Sample Tissue: Rat Spleen lysates
Antibody Dilution: 1 μ g/ml