

Product datasheet for **TA363303**

Hsd3b2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle terminal region of mouse HSD3B6
Specificity:	Expected reactivity: Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	42 kDa
Gene Name:	hydroxy-delta-5-steroid dehydrogenase, 3 beta- and steroid delta-isomerase 6
Database Link:	NP_058961.4 Entrez Gene 15493 Mouse O35469
Background:	enzyme responsible for the production of progesterone as well as the precursors of androgens, estrogens, glucocorticoids and mineralcorticoids [RGD, Feb 2006]
Synonyms:	D19210; Hsd3b1



[View online »](#)

Product images:



Host: Rabbit
Target Name: Hsd3b6
Sample Type: Mouse Testis Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: HSD3B6
Sample Tissue: Mouse Testis lysates
Antibody Dilution: 1 μ g/ml