

Product datasheet for **TA363326**

Anks6 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C terminal region of mouse ANKS6
Specificity:	Expected reactivity: Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	93 kDa
Gene Name:	ankyrin repeat and sterile alpha motif domain containing 6
Database Link:	NP_001019307.1 Entrez Gene 75691 Mouse Q6GQX6
Background:	Required for renal function.
Synonyms:	ANKRD14; DKFZp686D24121; DKFZp781I0117; FLJ36928; MGC70366; OTTHUMG00000020346; SamCystin; SAMD6



[View online »](#)

Product images:



Host: Rabbit
Target Name: Anks6
Sample Type: Mouse Liver Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: ANKS6
Sample Tissue: Mouse Liver lysates
Antibody Dilution: 1 μ g/ml